

## Product Datasheet

### Phytase Microplate Assay Kit (orb707341)

**Description** Detection and Quantification of Phytase Activity.

**Reactivity** All

**Range** 0.01 mmol/L - 1 mmol/L

**Tested Applications** FA

**Storage** Shipped and store at 4°C for 6 months.


**Note** For research use only

**Application notes** Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 660 nm


**Sample Types** Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

**CAS Number** 9001-89-2

**Expiration Date** Please enquire.



Representative Standard Curve



Positive Control reaction in 96-well pla...