

Product Datasheet

Sulfate Microplate Assay Kit (orb707342)

Description	Detection and Quantification of Sulfate Content.
Reactivity	All
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Other biological fluids, Tissue extracts, Cell lysate, Cell culture media, Platform: Microplate reader 420 nm
Sample Types	Urine, Serum, Plasma, Other biological fluids, Tissue extracts, Cell lysate, Cell culture media
CAS Number	14808-79-8
Expiration Date	Please enquire.