


Product Datasheet

Vitamin E Microplate Assay Kit (orb707345)

Description	Detection and Quantification of Vitamin E Content.
Reactivity	All
Range	50 µmol/L - 2000 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Other biological fluids, Platform: Microplate reader 530 nm
Sample Types	Serum, Plasma, Tissue extracts, Other biological fluids
CAS Number	10191-41-0
Expiration Date	Please enquire.



Representative Standard Curve