


Product Datasheet

3-alpha Hydroxysteroid Dehydrogenase Microplate Assay Kit (orb707346)

Description	Detection and Quantification of 3-alpha Hydroxysteroid Dehydrogenase Activity.
Reactivity	All
Range	4 µmol/L - 400 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Other biological fluids, Platform: Microplate reader 450 nm
Sample Types	Urine, Serum, Plasma, Other biological fluids
CAS Number	9028-56-2
Expiration Date	Please enquire.



Representative Standard Curve