


Product Datasheet

Lactate Microplate Assay Kit (orb707348)

Description	Detection and Quantification of Lactate Content.
Reactivity	All
Range	0.1 mmol/L - 100 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Cell culture media, Other biological fluids, Platform: Microplate reader 450 nm
Sample Types	Serum, Plasma, Cell culture media, Other biological fluids
CAS Number	10326-41-7
Expiration Date	Please enquire.



Representative Standard Curve