


Product Datasheet

Creatine Microplate Assay Kit (orb707351)

Description	Detection and Quantification of Creatine Content.
Reactivity	All
Range	0.05 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Other biological fluids, Platform: Microplate reader 546 nm
Sample Types	Urine, Serum, Plasma, Other biological fluids
CAS Number	6020-87-7
Expiration Date	Please enquire.



Representative Standard Curve