

Product Datasheet

Nitrite Microplate Assay Kit (orb707355)

Description Detection and Quantification of Nitrite Content.

Reactivity All

Range 0.01 nmol/ml - 1 nmol/ml

Tested Applications FA


Storage Shipped and store at 4°C for 6 months.

Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 547 nm

Sample Types Serum, Plasma, Urine, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

Expiration Date Please enquire.



Representative Standard Curve