


Product Datasheet

Albumin (Plant) Microplate Assay Kit (orb707358)

Description	Detection and Quantification of Albumin (Plant) Content.
Reactivity	All
Range	2 µg/ml - 200 µg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Powder, Platform: Microplate reader 595 nm
Sample Types	Tissue extracts, Powder
Expiration Date	Please enquire.



Representative
Standard
Curve