


Product Datasheet

alpha-L-Fucosidase Microplate Assay Kit (orb707364)

| | |
|----------------------------|--|
| Description | Detection and Quantification of alpha-L-Fucosidase Activity. |
| Reactivity | All |
| Range | 3 µmol/L - 300 µmol/L |
| Tested Applications | FA |
| Storage | Shipped and store at 4°C for 6 months. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 405 nm |
| Sample Types | Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids |
| CAS Number | 9037-65-4 |
| Expiration Date | Please enquire. |



Representative Standard Curve