


Product Datasheet

Selenium Microplate Assay Kit (orb707369)

Description	Detection and Quantification of Selenium Content.
Reactivity	All
Range	0.03 mmol/L - 0.3 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 450 nm
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Expiration Date	Please enquire.



Representative Standard Curve