


## Product Datasheet

### Caspase-3 Microplate Assay Kit (orb707372)

<b>Description</b>	Detection and Quantification of Caspase-3 Activity.
<b>Reactivity</b>	All
<b>Range</b>	0.01 µmol/ml - 5 µmol/ml
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Other biological fluids, Platform: Microplate reader 405 nm
<b>Sample Types</b>	Tissue extracts, Cell lysate, Other biological fluids
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve