


Product Datasheet

Blood Urea Nitrogen Microplate Assay Kit (orb707378)

Description	Detection and Quantification of Blood Urea Nitrogen Content.
Reactivity	All
Range	0.05 mmol/L - 1 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma Other biological fluids, Platform: Microplate reader 525 nm
Sample Types	Serum, Plasma Other biological fluids
Expiration Date	Please enquire.



Representative Standard Curve