


## Product Datasheet

### Creatinine Microplate Assay Kit (orb707380)

<b>Description</b>	Detection and Quantification of Creatinine Content.
<b>Reactivity</b>	All
<b>Range</b>	0.05 mmol/L - 5 mmol/L
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Other biological fluids, Platform: Microplate reader 546 nm
<b>Sample Types</b>	Urine, Serum, Plasma, Other biological fluids
<b>CAS Number</b>	60-27-5
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve