

## Product Datasheet

### Pectin Microplate Assay Kit (orb707382)

**Description** Detection and Quantification of Pectin Content.

**Reactivity** All

**Range** 0.3 mmol/L - 3 mmol/L

**Tested Applications** FA

**Storage** Shipped and store at 4°C for 6 months.


**Note** For research use only

**Application notes** Detection Method: Colorimetric, Sample Types: Tissue extracts, Other biological fluids, Platform: Microplate reader 540 nm

**Sample Types** Tissue extracts, Other biological fluids

**CAS Number** 9000-69-5

**Expiration Date** Please enquire.



Representative Standard Curve