


Product Datasheet

Phytic Acid Microplate Assay Kit (orb707384)

Description	Detection and Quantification of Phytic Acid Content.
Reactivity	All
Range	0.04 mmol/L - 4 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 660 nm
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	83-86-3
Expiration Date	Please enquire.



Representative Standard Curve