


Product Datasheet

Neutral Phosphatase Microplate Assay Kit (orb707385)

| | |
|----------------------------|---|
| Description | Detection and Quantification of Acid Phosphatase Activity. |
| Reactivity | All |
| Range | 0.05 mmol/L - 5 mmol/L |
| Tested Applications | FA |
| Storage | Shipped and store at 4°C for 6 months. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 510 nm |
| Sample Types | Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids |
| CAS Number | 9001-77-8 |
| Expiration Date | Please enquire. |



Representative Standard Curve