


Product Datasheet

Carboxylesterase Microplate Assay Kit (orb707386)

Description	Detection and Quantification of Carboxylesterase Activity.
Reactivity	All
Range	2 µmol/L - 200 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 600 nm
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Expiration Date	Please enquire.



Representative Standard Curve