



## Product Datasheet

# Glucose Dehydrogenase Microplate Assay Kit (orb707388)

<b>Description</b>	Detection and Quantification of Glucose Dehydrogenase Activity.
<b>Reactivity</b>	All
<b>Range</b>	4 µmol/L - 400 µmol/L
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 340 nm
<b>Sample Types</b>	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
<b>CAS Number</b>	9028-53-9
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve



Positive Control reaction in 96-well pla...