


Product Datasheet

Flavone Microplate Assay Kit (orb707390)

Description	Detection and Quantification of Flavone Content.
Reactivity	All
Range	0.02 mmol/L - 2 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Other biological fluids, Platform: Microplate reader 420 nm
Sample Types	Tissue extracts, Other biological fluids
CAS Number	525-82-6
Expiration Date	Please enquire.



Representative Standard Curve