


Product Datasheet

Flavonoid Microplate Assay Kit (orb707391)

Description	Detection and Quantification of Flavonoid Content.
Reactivity	All
Range	0.05 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Other biological fluids, Platform: Microplate reader 420 nm
Sample Types	Tissue extracts, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve