


Product Datasheet

Proanthocyanidin Microplate Assay Kit (orb707392)

Description	Detection and Quantification of Proanthocyanidin Content.
Reactivity	All
Range	0.1 mmol/L - 10 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Other biological fluids, Platform: Microplate reader 500 nm
Sample Types	Tissue extracts, Other biological fluids
CAS Number	4852-22-6
Expiration Date	Please enquire.



Representative Standard Curve