


## Product Datasheet

### Ethanol Microplate Assay Kit (orb707393)

<b>Description</b>	Detection and Quantification of Ethanol Content.
<b>Reactivity</b>	All
<b>Range</b>	2 µg/L - 200 µg/L
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine, Saliva, Other biological fluids, Platform: Microplate reader 450 nm
<b>Sample Types</b>	Serum, Plasma, Urine, Saliva, Other biological fluids
<b>CAS Number</b>	64-17-5
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve