

Product Datasheet

Lipoxygenase Microplate Assay Kit (orb707395)

Description Detection and Quantification of Lipoxygenase Activity.

Reactivity All

Range 0.03 mmol/L - 3 mmol/L

Tested Applications FA


Storage Shipped and store at 4°C for 6 months.

Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 470 nm

Sample Types Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

Expiration Date Please enquire.



Representative Standard Curve