


Product Datasheet

Total Iron-binding Capacity Microplate Assay Kit (orb707397)

| | |
|----------------------------|---|
| Description | Detection and Quantification of Total Iron-binding Capacity. |
| Reactivity | All |
| Range | 0.1 mmol/L - 10 mmol/L |
| Tested Applications | FA |
| Storage | Shipped and store at 4°C for 6 months. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric, Sample Types: Serum, Plasma, Platform: Microplate reader 562 nm |
| Sample Types | Serum, Plasma |
| Expiration Date | Please enquire. |



Representative Standard Curve