


Product Datasheet

Fructosamine Microplate Assay Kit (orb707398)

Description	Detection and Quantification of Fructosamine Content.
Reactivity	All
Range	0.01 mmol/L - 2 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Other biological fluids, Platform: Microplate reader 540 nm
Sample Types	Serum, Plasma, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve