


Product Datasheet

Tryptophan Microplate Assay Kit (orb707399)

Description	Detection and Quantification of Tryptophan Content.
Reactivity	All
Range	0.05 mmol/L - 10 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 600 nm
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	73-22-3
Expiration Date	Please enquire.



Representative
Standard
Curve