

Product Datasheet

Taurine Microplate Assay Kit (orb707400)

Description Detection and Quantification of Taurine Content.

Reactivity All

Range 0.1 mmol/L - 10 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 400 nm

Sample Types Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 107-35-7

Expiration Date Please enquire.



Representative
Standard
Curve