

Product Datasheet

Tuberin (phospho-Ser1387) antibody (orb7132)

Description	Tuberin (phospho-Ser1387) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Tuberin around the phosphorylation site of Ser1387 (SS(p-S)SP)
Target	Phospho-Tuberin (Ser1387)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	200 kDa
Uniprot ID	P49815
Expiration Date	12 months from date of receipt.