

Product Datasheet

Goat Anti-Rabbit IgG Secondary Antibody Biotin Conjugated (orb713463)

Description	Goat Anti-Rabbit IgG Secondary Antibody Biotin Conjugated
Species/Host	Goat
Conjugation	Biotin
Storage	Store at -20°C/1 year
Note	For research use only
Application notes	Working with enzyme-conjugated streptavidin ELISA: 1:20,000 - 1:200,000 Western Blot: 1:20,000 - 1:200,000 IHC/ICC: 1:500 - 1:5,000 Flow Cytometry: 1:200 - 1:1,000 IF 1:200 - 1:1,000
Clonality	Polyclonal
Dilution Range	Working with enzyme-conjugated streptavidin ELISA: 1:20,000 - 1:200,000 Western Blot: 1:20,000 - 1:200,000 IHC/ICC: 1:500 - 1:5,000 Flow Cytometry: 1:200 - 1:1,000 IF 1:200 - 1:1,000
Expiration Date	12 months from date of receipt.