

Product Datasheet

Goat Anti-Rabbit IgG Secondary Antibody Cy3 Conjugated (orb713464)

Description	Goat Anti-Rabbit IgG Secondary Antibody Cy3 Conjugated
Species/Host	Goat
Conjugation	Cy3
Storage	Store at -20°C/1 year
Note	For research use only
Application notes	Most Applications: 1:100 - 1:800
Clonality	Polyclonal
Dilution Range	Most Applications: 1:100 - 1:800
Expiration Date	12 months from date of receipt.