

## Product Datasheet

# Goat Anti-Rabbit IgG Secondary Antibody PE Conjugated (orb713466)

<b>Description</b>	Goat Anti-Rabbit IgG Secondary Antibody PE Conjugated
<b>Species/Host</b>	Goat
<b>Conjugation</b>	RPE
<b>Storage</b>	Store at -20°C/1 year
<b>Note</b>	For research use only
<b>Application notes</b>	Most Applications: 1:100 - 1:200
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	Most Applications: 1:100 - 1:200
<b>Expiration Date</b>	12 months from date of receipt.