

## Product Datasheet

### Batroxobin Antibody Blocking Peptide (orb72378)

<b>Description</b>	Batroxobin peptide.
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb11591 (Batroxobin antibody)
<b>Uniprot ID</b>	<b>P04971</b>
<b>Expiration Date</b>	12 months from date of receipt.