

## Product Datasheet

### Caspase 13 peptide (orb72426)

<b>Description</b>	Caspase 13 peptide.
<b>Target</b>	CASP13
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb10236 (Caspase 13/4 antibody)
<b>Uniprot ID</b>	<a href="#">O75601</a>
<b>Expiration Date</b>	12 months from date of receipt.