

## Product Datasheet

### Estrogen Receptor beta peptide (orb72657)

<b>Description</b>	Estrogen Receptor beta antibody
<b>Target</b>	ESR2
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb10616 (Estrogen Receptor beta antibody)
<b>Uniprot ID</b>	<a href="#">Q92731</a>
<b>Expiration Date</b>	12 months from date of receipt.