

## Product Datasheet

### Smad2 (phospho-Ser465) peptide (orb73214)

<b>Description</b>	Smad2 (phospho-Ser465) peptide.
<b>Target</b>	SMAD2
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q15796</a>
<b>Expiration Date</b>	12 months from date of receipt.