

Product Datasheet

Thymosin alpha-1 peptide (orb73269)

Description	Thymosin alpha-1 peptide.
Target	PTMA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb11478 (Thymosin alpha 1 antibody)
Protein Sequence	Ac-Ser-Asp-Ala-Ala-Val-Asp-Thr-Ser-Ser-Glu-Ile-Thr-Thr-Lys-Asp-Leu-Lys-Glu-Lys-Lys-Glu-Val-Val-Glu-Glu-Ala-Glu-Asn
MW	3108 kDa
Uniprot ID	P06454
Expiration Date	12 months from date of receipt.