



Product Datasheet

ATPase beta3(Na⁺/K⁺) antibody (orb73384)

| | |
|----------------------------|---|
| Description | ATPase beta3(Na ⁺ /K ⁺) antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P54709 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of A:Hela whole ce...



Immunohistochemical analysis of formalin...