



## Product Datasheet

### Cortisol Binding Globulin antibody (orb73445)


|                            |   |
|----------------------------|---|
| <b>Description</b>         | Cortisol Binding Globulin antibody  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IF, WB  |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Recombinant   |
| <b>Uniprot ID</b>          | <b>P08185</b>   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Western Blot analysis of A:A431 whole ce...



Merged with DNA probe



Immunofluorescence analysis of formalin-

...