


Product Datasheet

Plexin-D1 antibody (orb73475)

Description	Plexin-D1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y4D7
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:HepG2 whole c...



Confocal immunofluorescence analysis of ...