

## Product Datasheet

### Human NOTCH2 Protein (orb735046)

<b>Description</b>	Human NOTCH2, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Gln 530 (Ac...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	NOTCH2
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 80.5 kDa. The protein migrates as 90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_077719.2
<b>Purity</b>	95%
<b>Source</b>	Human NOTCH2, Fc Tag (orb735046) is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Gln 530 (Accession # Q04721-1).
<b>MW</b>	80.5 kDa
<b>Uniprot ID</b>	<a href="#">Q04721</a>
<b>NCBI</b>	<a href="#">NP_077719.2</a>
<b>Expiration Date</b>	6 months from date of receipt.