



Product Datasheet

NKG2-A (CD159a) antibody (orb73560)

Description	NKG2-A (CD159a) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P26715
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431, B:HepG2...



Confocal immunofluorescence analysis of ...