



Product Datasheet

SH3 domain-binding protein 1 antibody (orb73594)

Description	SH3 domain-binding protein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y3L3
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Hela, B:Jurka...



Costained with Hoechst 33342

Confocal immunofluorescence analysis of ...