



Product Datasheet

PPM1K antibody (orb73604)

Description	PPM1K antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8N3J5
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293T cell line ...



Immunohistochemical analysis of formalin...