



Product Datasheet

RGS17 antibody (orb73627)


Description	RGS17 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9UGC6
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:HepG2 whole c...



Costained with Hoechst 33342



Confocal immunofluorescence analysis of ...