



Product Datasheet

CD55 antibody (orb73639)

Description	CD55 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08174
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western Blot analysis of A:THP-1 whole c...