

Product Datasheet

CD93 antibody (orb736702)

Description	Mouse monoclonal antibody to CD93
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Target	CD93
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	68560
Uniprot ID	Q9NPY3
Entrez	22918
Expiration Date	12 months from date of receipt.