

Product Datasheet

PTGS2 antibody (orb736751)

Description	Rabbit monoclonal antibody to PTGS2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Target	PTGS2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	68996
Uniprot ID	P35354
Entrez	5743
Dilution Range	IHC: Ready-to-use
Expiration Date	12 months from date of receipt.