

Product Datasheet

HLA-DRB1 antibody (orb736790)

Description	Mouse monoclonal antibody to HLA-DRB1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Target	HLA-DRB1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
MW	29966
Uniprot ID	P01911
Entrez	3123
Dilution Range	IHC: Ready-to-use
Expiration Date	12 months from date of receipt.