

Product Datasheet

RBFOX3 antibody (orb736821)

Description	Rabbit monoclonal antibody to RBFOX3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Target	RBFOX3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	33873
Uniprot ID	A6NFN3
Entrez	146713
Dilution Range	IHC: Ready-to-use
Expiration Date	12 months from date of receipt.