

Product Datasheet

TAG-72 antibody (orb736862)

Description	Mouse monoclonal antibody to TAG-72
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	-
Dilution Range	IHC: Ready-to-use
Expiration Date	12 months from date of receipt.